

PiroGon Device API

1.0.0

OAS 3.1

[/openapi.json](#)

Local-based PiroGon device data management system

Authorize



authentication



POST

/api/v1/auth/login Login



User login with email and password

Parameters

Try it out

Reset

No parameters

Request body required

application/json



Example Value | Schema

```
{
  "email": "superadmin@pirogon.com",
  "password": "SuperSecret123!"
}
```

Responses

Curl

```
curl -X 'POST' \
  'http://127.0.0.1:8000/api/v1/auth/login' \
  -H 'accept: application/json' \
  -H 'Content-Type: application/json' \
```

```
-d '{
  "email": "superadmin@pirogon.com",
  "password": "SuperSecret123!"
}'
```

Request URL

```
http://127.0.0.1:8000/api/v1/auth/login
```

Server response

Code	Details
------	---------

200

Response body

```
{
  "access_token": "eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1YzhINS1iYzI3LTRmZmYtYWMyNy1jOGY1Mzh1MTQ5MzYiLCJlbWFPbCI6InN1cGVyYWRtaW5AcGlyb2d2b20iLCJyb2x1IjoieWVhbnRtaW51LCJleHAiOjE3NTg5MTc2ODQ5Im1hdCI6MTc1ODkxNTg4NH0uYjFkY2F7XxqgwPD21P7i9K98oq7DSErhMKFM20IbpTA",
  "token_type": "bearer",
  "expires_in": 1800,
  "user": {
    "id": "0435c8e5-bc27-4fff-ac67-c8f538e14936",
    "email": "superadmin@pirogon.com",
    "full_name": "System Administrator",
    "role": "admin"
  }
}
```

 [Download](#)

Response headers

```
access-control-allow-credentials: true
content-length: 447
content-type: application/json
date: Fri, 26 Sep 2025 19:44:43 GMT
server: uvicorn
```

Responses

Code	Description	Links
200	Successful Response	No links
	Media type application/json Controls Accept header. Example Value Schema<pre>{ "access_token": "string", "token_type": "bearer", "expires_in": 0, "user": {} }</pre>	
422	Validation Error	No links
	Media type application/json	

Code

Description

Links

Example Value

Schema

```
{  "detail": [    {      "loc": [        "string",        0      ],      "msg": "string",      "type": "string"    }  ]}
```

POST

/api/v1/auth/kiosk

Kiosk Login

^

Kiosk login with PIN

Parameters

Try it out

No parameters

Request body required

application/json

Example Value

Schema

```
{  "pin": "string"}
```

Responses

Code

Description

Links

200

Successful Response

No links

Media type

application/json

Controls Accept header.

Example Value

Schema

Code	Description	Links
	<pre>{ "access_token": "string", "token_type": "bearer", "expires_in": 0, "user": {}}</pre>	
422	Validation Error Media type application/json Example Value Schema <pre>{ "detail": [{ "loc": ["string", 0], "msg": "string", "type": "string" }]}</pre>	No links

GET

/api/v1/auth/me

Get Current User Info



Get current user information

Parameters

Try it out

No parameters

Responses

Curl

```
curl -X 'GET' \
'http://127.0.0.1:8000/api/v1/auth/me' \
-H 'accept: application/json' \
-H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1YzhINS1iYzI3IHRm'
```



Request URL

http://127.0.0.1:8000/api/v1/auth/me

Server response

Code	Details
200	Response body <pre>{ "id": "0435c8e5-bc27-4fff-ac67-c8f538e14936", "email": "superadmin@pirogon.com", "full_name": "System Administrator", "role": "admin", "is_active": true, "last_login": "2025-09-26T19:03:28.190620" }</pre>  Download Response headers <pre>content-length: 187 content-type: application/json date: Fri, 26 Sep 2025 19:04:14 GMT server: uvicorn</pre>

Responses

Code	Description	Links
200	Successful Response Media type application/json ▼ Controls Accept header. Example Value Schema "string"	No links

POST

/api/v1/auth/create-user

Create User

 

Create a new user - Admin only

Parameters

Try it out

Reset

No parameters

Request body required

application/json

▼

Example Value | Schema

```
{
  "email": "operator1@pirogon.com",
  "username": "operator1",
  "full_name": "First Operator",
  "password": "OperatorPass123!",
  "role": "operator"
}
```

Responses

Curl

```
curl -X 'POST' \
  'http://127.0.0.1:8000/api/v1/auth/create-user' \
  -H 'accept: application/json' \
  -H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1Yzh1NS1iYzI3LTRm' \
  -H 'Content-Type: application/json' \
  -d '{
    "email": "operator1@pirogon.com",
    "username": "operator1",
    "full_name": "First Operator",
    "password": "OperatorPass123!",
    "role": "operator"
  }'
```

Request URL

```
http://127.0.0.1:8000/api/v1/auth/create-user
```

Server response

Code

Details

200

Response body

```
{
  "id": "32699633-9cfb-47aa-b23d-9634ddf75bac",
  "email": "operator1@pirogon.com",
  "username": "operator1",
  "full_name": "First Operator",
  "role": "operator",
  "is_active": true
}
```



Download

Response headers

```
access-control-allow-credentials: true
content-length: 164
content-type: application/json
date: Fri, 26 Sep 2025 19:06:16 GMT
server: uvicorn
```

Responses

Code	Description	Links
200	Successful Response	No links
Media type application/json Controls Accept header. Example Value Schema <pre>{}</pre>		
422	Validation Error	No links
Media type application/json Example Value Schema <pre>{ "detail": [{ "loc": ["string", 0], "msg": "string", "type": "string" }]}</pre>		

patients

^

GET

/api/v1/patients/

Get Patients

Get all patients with optional search and pagination

Parameters

Try it out

Name	Description
skip integer (query)	Default value : 0 0 minimum: 0

Code

Details

```
"patient_number": "P1004",
"name": "John4 Doe4",
"date_of_birth": "1997-07-21T00:00:00",
"sex": "male",
"medical_record_number": "MRN-42345",
"id": "1663e6f6-8d3b-4a50-8be2-d3d4241ebf02"
```



Download

Response headers

```
content-length: 1049
content-type: application/json
date: Fri,26 Sep 2025 19:45:09 GMT
server: uvicorn
```

Responses

Code	Description	Links
200	Successful Response	No link
Media type application/json Controls Accept header. Example Value Schema [{ "patient_number": "string", "name": "string", "date_of_birth": "2025-09-26T20:12:00.083Z", "sex": "male", "medical_record_number": "string", "id": "string", "created_at": "2025-09-26T20:12:00.083Z", "updated_at": "2025-09-26T20:12:00.083Z" }]		
422	Validation Error	No links
Media type application/json Example Value Schema { "detail": [{ "loc": ["string", 0], "msg": "string", "type": "string" }] }		

POST

/api/v1/patients/ Create Patient



Create a new patient

Parameters

Try it out

Reset

No parameters

Request body required

application/json

Example Value

Schema

```
{
  "patient_number": "P1001",
  "name": "John Doe",
  "date_of_birth": "1997-07-21",
  "sex": "male",
  "medical_record_number": "MRN-12345"
}
```

Responses

Curl

```
curl -X 'POST' \
  'http://127.0.0.1:8000/api/v1/patients/' \
  -H 'accept: application/json' \
  -H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1YzhINS1iYzI3LTRm
  -H 'Content-Type: application/json' \
  -d '{
    "patient_number": "P1001",
    "name": "John Doe",
    "date_of_birth": "1997-07-21",
    "sex": "male",
    "medical_record_number": "MRN-12345"
  }'
```



Request URL

http://127.0.0.1:8000/api/v1/patients/

Server response

Code	Details
201	Response body

Code

Details

```
{
  "patient_number": "P1001",
  "name": "John Doe",
  "date_of_birth": "1997-07-21T00:00:00",
  "sex": "male",
  "medical_record_number": "MRN-12345",
  "id": "2707a3aa-c798-42af-a94e-40e02ed18a1d",
  "created_at": "2025-09-26T19:29:20.938096",
  "updated_at": "2025-09-26T19:29:20.938101"
}
```

 [Download](#)

Response headers

```
access-control-allow-credentials: true
content-length: 259
content-type: application/json
date: Fri, 26 Sep 2025 19:29:20 GMT
server: uvicorn
```

Responses

Code	Description	Links
201	Successful Response	No links
Media type application/json Controls Accept header. Example Value Schema <pre>{ "patient_number": "string", "name": "string", "date_of_birth": "2025-09-26T20:12:00.086Z", "sex": "male", "medical_record_number": "string", "id": "string", "created_at": "2025-09-26T20:12:00.086Z", "updated_at": "2025-09-26T20:12:00.086Z" }</pre>		
422	Validation Error	No links
Media type application/json Example Value Schema <pre>{ "detail": [{ "loc": ["string", 0], "msg": "string", "type": "string" }] }</pre>		

Code

Description

Links

```
  ]
}
```

GET

/api/v1/patients/{patient_id} Get Patient



Get a specific patient by ID

Parameters

Try it out

Name

Description

patient_id * required

string

(path)

fdd96bb7-08b2-49fb-ad7e-a0f566de59a6

Responses

Curl

```
curl -X 'GET' \
'http://127.0.0.1:8000/api/v1/patients/fdd96bb7-08b2-49fb-ad7e-a0f566de59a6' \
-H 'accept: application/json' \
-H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1Yzh1NS1iYzI3IiwiaWF0Ij0iMTY5OTUyMjY0IiwiaXNjaWkiOiJ1b3RlciJ9.eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1Yzh1NS1iYzI3IiwiaWF0Ij0iMTY5OTUyMjY0IiwiaXNjaWkiOiJ1b3RlciJ9'
```

Request URL

http://127.0.0.1:8000/api/v1/patients/fdd96bb7-08b2-49fb-ad7e-a0f566de59a6

Server response

Code

Details

200

Response body

```
{
  "patient_number": "P1001",
  "name": "John Doe updated",
  "date_of_birth": "1985-07-21T00:00:00",
  "sex": "male",
  "medical_record_number": "MRN-12345",
  "id": "fdd96bb7-08b2-49fb-ad7e-a0f566de59a6",
  "created_at": "2025-09-26T19:08:37.514710",
  "updated_at": "2025-09-26T19:13:16.368051"
}
```



Download

Code

Details

Response headers

```
content-length: 267
content-type: application/json
date: Fri,26 Sep 2025 19:13:25 GMT
server: uvicorn
```

Responses

Code	Description	Links
------	-------------	-------

200	Successful Response	No links
-----	---------------------	----------

Media type

application/json

Controls Accept header.

Example Value | Schema

```
{
  "patient_number": "string",
  "name": "string",
  "date_of_birth": "2025-09-26T20:12:00.088Z",
  "sex": "male",
  "medical_record_number": "string",
  "id": "string",
  "created_at": "2025-09-26T20:12:00.088Z",
  "updated_at": "2025-09-26T20:12:00.088Z"
}
```

422	Validation Error	No links
-----	------------------	----------

Media type

application/json

Example Value | Schema

```
{
  "detail": [
    {
      "loc": [
        "string",
        0
      ],
      "msg": "string",
      "type": "string"
    }
  ]
}
```

PUT

/api/v1/patients/{patient_id} Update Patient



Update a patient

Parameters

Try it out

Reset

NameDescription

patient_id * required

string
(path)

2707a3aa-c798-42af-a94e-40e02ed18a1d

Request body required

application/json

Example Value Schema

```
{
  "patient_number": "P1001",
  "name": "John Doe updated",
  "date_of_birth": "1985-07-21T00:00:00",
  "sex": "male",
  "medical_record_number": "MRN-12345"
}
```

Responses

Curl

```
curl -X 'PUT' \
  'http://127.0.0.1:8000/api/v1/patients/2707a3aa-c798-42af-a94e-40e02ed18a1d' \
  -H 'accept: application/json' \
  -H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1YzhINS1iYzI3LTRm' \
  -H 'Content-Type: application/json' \
  -d '
{
  "patient_number": "P1001",
  "name": "John Doe updated",
  "date_of_birth": "1985-07-21T00:00:00",
  "sex": "male",
  "medical_record_number": "MRN-12345"
}'
```

Request URL

http://127.0.0.1:8000/api/v1/patients/2707a3aa-c798-42af-a94e-40e02ed18a1d

Server response

Code

Details

500

Undocumented

Error: Internal Server Error

Response body

Code

Details

Internal Server Error

 Download

Response headers

```
content-length: 21
content-type: text/plain; charset=utf-8
date: Fri,26 Sep 2025 19:45:47 GMT
server: uvicorn
```

Responses

Code	Description	Links
200	Successful Response	No links

Media type

application/json

Controls Accept header.

Example Value | Schema

```
{
  "patient_number": "string",
  "name": "string",
  "date_of_birth": "2025-09-26T20:12:00.092Z",
  "sex": "male",
  "medical_record_number": "string",
  "id": "string",
  "created_at": "2025-09-26T20:12:00.092Z",
  "updated_at": "2025-09-26T20:12:00.092Z"
}
```

422	Validation Error	No links
-----	------------------	----------

Media type

application/json

Example Value | Schema

```
{
  "detail": [
    {
      "loc": [
        "string",
        0
      ],
      "msg": "string",
      "type": "string"
    }
  ]
}
```

DELETE /api/v1/patients/{patient_id} Delete Patient



Parameters

Name

patient_id * required
string
(path)

```
curl -X 'DELETE' \
'http://127.0.0.1:8000/api/v1/patients/fdd96bb7-08b2-49fb-ad7e-a0f566de59a6' \
-H 'accept: application/json' \
-H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1Yzh1NS1iYzI3LTRm'
```

<http://127.0.0.1:8000/api/v1/patients/fdd96bb7-08b2-49fb-ad7e-a0f566de59a6>

Code

200

```
{
  "message": "Patient deleted successfully"
}
```

```
access-control-allow-credentials: true
content-length: 42
content-type: application/json
date: Fri, 26 Sep 2025 19:28:17 GMT
server: uvicorn
```

Code

Links

200

Successful Response

No links

Code	Description	Links
	Media type application/json Controls Accept header. Example Value Schema "string"	
422	Validation Error Media type application/json Example Value Schema { "detail": [{ "loc": ["string", 0], "msg": "string", "type": "string" }] }	No links

sessions

GET	/api/v1/sessions/	Get Sessions		
Get sessions with optional filtering				
Parameters				Try it out
Name	Description			
patient_id	2707a3aa-c798-42af-a94e-40e02ed18a1d			
string (string null) (query)				

Name	Description
------	-------------

device_id string (string null) (query)	28188113-7d98-4e99-b96e-5574d2baf69f
--	---

skip integer (query)	Default value : 0 0
----------------------------	------------------------------------

limit integer (query)	Default value : 100 100
-----------------------------	--

Responses

Curl

```
curl -X 'GET' \  
  'http://127.0.0.1:8000/api/v1/sessions/?patient_id=2707a3aa-c798-42af-a94e-40e02ed18a1d&device_  
  -H 'accept: application/json' \  
  -H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1YzhINS1iYzI3IHRm
```

Request URL

```
http://127.0.0.1:8000/api/v1/sessions/?patient_id=2707a3aa-c798-42af-a94e-  
40e02ed18a1d&device_id=28188113-7d98-4e99-b96e-5574d2baf69f&skip=0&limit=100
```

Server response

Code	Details
------	---------

200	Response body [{ "session_number": 1, "threshold_level": 10, "tolerance_level": 10, "stimulation_level": 10, "frequency_level": 20, "pulse_width_level": 200, "operator_name": "John Doe", "notes": "string", "id": "a0d09fc1-5386-4aac-ba64-fb4fd5e0dae5", "start_time": "2025-09-26T19:46:03.997468", "end_time": null, "status": "active", "patient_id": "2707a3aa-c798-42af-a94e-40e02ed18a1d", "device_id": "28188113-7d98-4e99-b96e-5574d2baf69f", "operator_type": "login", "test_passed": null, }]
-----	--

Code

Details



Download

```
"stimulation_status": null
}
]
```

Response headers

```
content-length: 478
content-type: application/json
date: Fri, 26 Sep 2025 19:47:25 GMT
server: uvicorn
```

Responses

Code	Description	Links
200	Successful Response	No links

Media type

application/json

Controls Accept header.

Example Value | Schema

```
[
{
  "session_number": 1,
  "threshold_level": 10,
  "tolerance_level": 10,
  "stimulation_level": 10,
  "frequency_level": 1,
  "pulse_width_level": 1,
  "operator_name": "string",
  "notes": "string",
  "id": "string",
  "start_time": "2025-09-26T20:12:00.098Z",
  "end_time": "2025-09-26T20:12:00.098Z",
  "status": "active",
  "patient_id": "string",
  "device_id": "string",
  "operator_type": "login",
  "test_passed": true,
  "stimulation_status": "string"
}
]
```

422	Validation Error	No links
-----	------------------	----------

Media type

application/json

Example Value | Schema

```
{
  "detail": [
    {
      "loc": [
        "string",
        0
      ],
      "msg": "string",
    }
  ]
}
```

Code

Description

Links

```
      "type": "string"
    }
  ]
}
```

POST

/api/v1/sessions/ Create Session



Create a new therapy session

Parameters

Try it out

Reset

No parameters

Request body required

application/json 

Example Value

Schema

```
{
  "session_number": 1,
  "threshold_level": 10,
  "tolerance_level": 25,
  "stimulation_level": 17.5,
  "frequency_level": 20,
  "pulse_width_level": 200,
  "operator_name": "John Doe",
  "notes": "string",
  "patient_id": "2707a3aa-c798-42af-a94e-40e02ed18a1d",
  "device_id": "28188113-7d98-4e99-b96e-5574d2baf69f",
  "operator_type": "login"
}
```

Responses

Curl

```
curl -X 'POST' \
  'http://127.0.0.1:8000/api/v1/sessions/' \
  -H 'accept: application/json' \
  -H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1Yzh1NS1iYzI3LTRm
  -H 'Content-Type: application/json' \
  -d '{
    "session_number": 1,
    "threshold_level": 10,
    "tolerance_level": 25,
    "stimulation_level": 17.5,
    "frequency_level": 20,
    "pulse_width_level": 200,
```

```
"operator_name": "John Doe",
"notes": "string",
"patient_id": "2707a3aa-c798-42af-a94e-40e02ed18a1d",
"device_id": "28188113-7d98-4e99-b96e-5574d2baf69f",
"operator_type": "login"
},'
```



Request URL

```
http://127.0.0.1:8000/api/v1/sessions/
```

Server response

Code Details

201

Response body

```
{
  "session_number": 1,
  "threshold_level": 10,
  "tolerance_level": 25,
  "stimulation_level": 17.5,
  "frequency_level": 20,
  "pulse_width_level": 200,
  "operator_name": "John Doe",
  "notes": "string",
  "id": "4cd4431a-64fb-437f-b8d2-b6ae12ca69d2",
  "start_time": "2025-09-26T20:05:51.181609",
  "end_time": null,
  "status": "active",
  "patient_id": "2707a3aa-c798-42af-a94e-40e02ed18a1d",
  "device_id": "28188113-7d98-4e99-b96e-5574d2baf69f",
  "operator_type": "login",
  "test_passed": null,
  "stimulation_status": null
}
```

Download

Response headers

```
access-control-allow-credentials: true
content-length: 476
content-type: application/json
date: Fri, 26 Sep 2025 20:05:50 GMT
server: uvicorn
```

Responses

Code Description Links

201 Successful Response No links

Media type

application/json ▼

Controls Accept header.

Example Value | Schema

```
{
  "session_number": 1,
  "threshold_level": 10,
  "tolerance_level": 10,
  "stimulation_level": 10,
```

Code	Description	Links
	<pre>"frequency_level": 1, "pulse_width_level": 1, "operator_name": "string", "notes": "string", "id": "string", "start_time": "2025-09-26T20:12:00.102Z", "end_time": "2025-09-26T20:12:00.102Z", "status": "active", "patient_id": "string", "device_id": "string", "operator_type": "login", "test_passed": true, "stimulation_status": "string" }</pre>	
422	Validation Error	No links
Media type		
application/json		
Example Value Schema		
<pre>{ "detail": [{ "loc": ["string", 0], "msg": "string", "type": "string" }] }</pre>		

GET

/api/v1/sessions/{session_id}

Get Session



Get a specific session

Parameters

Try it out

Name	Description
session_id * required	
string (path)	a0d09fc1-5386-4aac-ba64-fb4fd5e0dae5

Responses

Curl

```
curl -X 'GET' \
'http://127.0.0.1:8000/api/v1/sessions/a0d09fc1-5386-4aac-ba64-fb4fd5e0dae5' \
-H 'accept: application/json' \
-H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1Yzh1NS1iYzI3IHRm
```

Request URL

```
http://127.0.0.1:8000/api/v1/sessions/a0d09fc1-5386-4aac-ba64-fb4fd5e0dae5
```

Server response

Code	Details
200	Response body <pre>{ "session_number": 1, "threshold_level": 10, "tolerance_level": 10, "stimulation_level": 10, "frequency_level": 20, "pulse_width_level": 200, "operator_name": "John Doe", "notes": "string", "id": "a0d09fc1-5386-4aac-ba64-fb4fd5e0dae5", "start_time": "2025-09-26T19:46:03.997468", "end_time": "2025-09-26T19:56:33.526612", "status": "completed", "patient_id": "2707a3aa-c798-42af-a94e-40e02ed18a1d", "device_id": "28188113-7d98-4e99-b96e-5574d2baf69f", "operator_type": "login", "test_passed": true, "stimulation_status": "active" }</pre>  Download Response headers <pre>content-length: 507 content-type: application/json date: Fri,26 Sep 2025 19:56:44 GMT server: uvicorn</pre>

Responses

Code	Description	Links
200	Successful Response Media type application/json Controls Accept header. Example Value Schema <pre>{ "session_number": 1, "threshold_level": 10,</pre>	No links

Code	Description	Links
	<pre>"tolerance_level": 10, "stimulation_level": 10, "frequency_level": 1, "pulse_width_level": 1, "operator_name": "string", "notes": "string", "id": "string", "start_time": "2025-09-26T20:12:00.105Z", "end_time": "2025-09-26T20:12:00.105Z", "status": "active", "patient_id": "string", "device_id": "string", "operator_type": "login", "test_passed": true, "stimulation_status": "string" }</pre>	
422	Validation Error	No links
Media type		
application/json		
Example Value Schema		
<pre>{ "detail": [{ "loc": ["string", 0], "msg": "string", "type": "string" }] }</pre>		

PUT

/api/v1/sessions/{session_id}/complete

Complete Session



Complete a therapy session

Parameters

Try it out

Reset

Name	Description
session_id * required	
string	a0d09fc1-5386-4aac-ba64-fb4fd5e0dae5
(path)	

Request body required

application/x-www-form-urlencoded

test_passed * required

boolean

stimulation_status * required

string

notes

string | (string | null)

Responses

Curl

```
curl -X 'PUT' \
  'http://127.0.0.1:8000/api/v1/sessions/a0d09fc1-5386-4aac-ba64-fb4fd5e0dae5/complete' \
  -H 'accept: application/json' \
  -H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1Yzh1NS1iYzI3LTRm' \
  -H 'Content-Type: application/x-www-form-urlencoded' \
  -d 'test_passed=true&stimulation_status=active&notes=string'
```

Request URL

http://127.0.0.1:8000/api/v1/sessions/a0d09fc1-5386-4aac-ba64-fb4fd5e0dae5/complete

Server response

Code	Details
200	Response body <pre>{ "message": "Session completed successfully" }</pre> Download Response headers <pre>access-control-allow-credentials: true content-length: 44 content-type: application/json date: Fri, 26 Sep 2025 19:56:33 GMT server: uvicorn</pre>

Responses

Code	Description	Links
200	Successful Response	No links

Code	Description	Links
	Media type application/json Controls Accept header. Example Value Schema "string"	
422	Validation Error	No links
	Media type application/json Example Value Schema <pre>{ "detail": [{ "loc": ["string", 0], "msg": "string", "type": "string" }]}</pre>	

POST

/api/v1/sessions/{session_id}/upload/{signal_type}

Upload Signal Data



Upload signal data for a session

Parameters

Try it out

Name	Description
session_id * required string (path)	session_id
signal_type * required string (path)	signal_type

Request body required

multipart/form-data

▼

file * required
string(\$binary)

Responses

Code	Description	Links
200	Successful Response Media type application/json Controls Accept header. Example Value Schema "string"	No links
422	Validation Error Media type application/json Example Value Schema { "detail": [{ "loc": ["string", 0], "msg": "string", "type": "string" }] }	No links

GET

/api/v1/sessions/{session_id}/export

Export Session



Export session data as ZIP file

Parameters

Try it out

Name	Description
session_id * required string (path)	session_id

Responses

Code	Description	Links
200	Successful Response Media type application/json ▾ Controls Accept header. Example Value Schema "string"	No links

Code	Description	Links
422	Validation Error	No links

Media type

application/json

Example Value

Schema

```
{  "detail": [    {      "loc": [        "string",        0      ],      "msg": "string",      "type": "string"    }  ]}
```

devices



GET

/api/v1/devices/ Get Devices

^

Get all registered devices

Parameters

Try it out

No parameters

Responses

Curl

```
curl -X 'GET' \
'http://127.0.0.1:8000/api/v1/devices/' \
-H 'accept: application/json' \
-H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1Yzh1NS1iYzI3IHRm
```

Request URL

```
http://127.0.0.1:8000/api/v1/devices/
```

Server response

Code

Details

200

Response body

```
[
  {
    "id": "28188113-7d98-4e99-b96e-5574d2baf69f",
    "serial_number": "CM5-001",
    "name": "Medical Device CM5",
    "model": null,
    "firmware_version": "1.0.0",
    "software_version": "1.0.0",
    "status": "active",
    "last_calibration": null,
    "next_maintenance": null,
    "registered_at": "2025-09-26T18:31:11.614161"
  }
]
```



Download

Response headers

```
content-length: 279
content-type: application/json
date: Fri,26 Sep 2025 19:49:19 GMT
server: uvicorn
```

Responses

Code	Description	Links
200	Successful Response	No links

Media type

application/json

Controls Accept header.

Example Value | Schema

```
[
  {
    "id": "string",
    "serial_number": "string",
    "name": "string",
    "model": "string",
    "firmware_version": "string",
    "software_version": "string",
    "status": "active",
    "last_calibration": "2025-09-26T20:12:00.116Z",
    "next_maintenance": "2025-09-26T20:12:00.116Z",
    "registered_at": "2025-09-26T20:12:00.116Z"
  }
]
```

GET

/api/v1/devices/current Get Current Device



Get current device information

Parameters

Try it out

No parameters

Responses

Curl

```
curl -X 'GET' \
  'http://127.0.0.1:8000/api/v1/devices/current' \
  -H 'accept: application/json' \
  -H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1YzhINS1iYzI3IHRm
```

Request URL

```
http://127.0.0.1:8000/api/v1/devices/current
```

Server response

Code	Details
------	---------

200

Response body

```
{
  "id": "28188113-7d98-4e99-b96e-5574d2baf69f",
  "serial_number": "CM5-001",
  "name": "Medical Device CM5",
  "model": null,
  "firmware_version": "1.0.0",
  "software_version": "1.0.0",
  "status": "active",
  "last_calibration": null,
  "next_maintenance": null,
  "registered_at": "2025-09-26T18:31:11.614161"
}
```

Download

Response headers

```
content-length: 277
content-type: application/json
date: Fri,26 Sep 2025 19:58:03 GMT
server: uvicorn
```

Responses

Code	Description	Links
------	-------------	-------

200

Successful Response

No links

Media type

application/json

Controls Accept header.

Example Value | Schema

Code

Description

Links

```
{
  "id": "string",
  "serial_number": "string",
  "name": "string",
  "model": "string",
  "firmware_version": "string",
  "software_version": "string",
  "status": "active",
  "last_calibration": "2025-09-26T20:12:00.125Z",
  "next_maintenance": "2025-09-26T20:12:00.125Z",
  "registered_at": "2025-09-26T20:12:00.125Z"
}
```

PUT

/api/v1/devices/current/status

Update Device Status



Update current device status

Parameters

Try it out

Name	Description
status * required string (query)	active

Responses

Curl

```
curl -X 'PUT' \
  'http://127.0.0.1:8000/api/v1/devices/current/status?status=active' \
  -H 'accept: application/json' \
  -H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1Yzh1NS1iYzI3IHRm
```



Request URL

```
http://127.0.0.1:8000/api/v1/devices/current/status?status=active
```

Server response

Code	Details
200	Response body

Code

Details

```
{
  "message": "Device status updated successfully"
}
```

 [Download](#)

Response headers

```
access-control-allow-credentials: true
content-length: 48
content-type: application/json
date: Fri, 26 Sep 2025 19:57:48 GMT
server: uvicorn
```

Responses

Code	Description	Links
200	Successful Response	No links
Media type application/json ▾ Controls Accept header. Example Value Schema "string"		
422	Validation Error	No links
Media type application/json ▾ Example Value Schema { "detail": [{ "loc": ["string", 0], "msg": "string", "type": "string" }] }		

system

^

GET	/api/v1/system/health	Health Check	^
-----	-----------------------	--------------	---

System health check endpoint

Parameters

Try it out

No parameters

Responses

Curl

```
curl -X 'GET' \
  'http://127.0.0.1:8000/api/v1/system/health' \
  -H 'accept: application/json'
```

Request URL

```
http://127.0.0.1:8000/api/v1/system/health
```

Server response

Code	Details
200	Response body <pre>{ "status": "healthy", "timestamp": "2025-09-26T20:02:49.907312", "version": "1.0.0" }</pre> Response headers <pre>content-length: 79 content-type: application/json date: Fri, 26 Sep 2025 20:02:49 GMT server: uvicorn</pre>

Responses

Code	Description	Links
200	Successful Response Media type application/json Controls Accept header. Example Value Schema	No links

Code

Description

Links

"string"

GET

/api/v1/system/stats

Get System Stats



Get system statistics

Parameters

Try it out

No parameters

Responses

Curl

```
curl -X 'GET' \
  'http://127.0.0.1:8000/api/v1/system/stats' \
  -H 'accept: application/json' \
  -H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1Yzh1NS1iYzI3IHRm
```

Request URL

```
http://127.0.0.1:8000/api/v1/system/stats
```

Server response

Code

Details

200

Response body

```
{
  "patients": 4,
  "sessions": 2,
  "active_sessions": 1,
  "completed_sessions": 1
}
```



Download

Response headers

```
content-length: 70
content-type: application/json
date: Fri, 26 Sep 2025 20:06:48 GMT
server: uvicorn
```

Responses

Code	Description	Links
200	Successful Response	No links
Media type application/json Controls Accept header. Example Value Schema "string"		

POST

/api/v1/system/backup

Create Backup



Create system backup

Parameters

Try it out

No parameters

Responses

Curl

```
curl -X 'POST' \
  'http://127.0.0.1:8000/api/v1/system/backup' \
  -H 'accept: application/json' \
  -H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1Yzh1NS1iYzI3LTRm' \
  -d ''
```



Request URL

http://127.0.0.1:8000/api/v1/system/backup

Server response

Code	Details
200	Response body <pre>{ "message": "Backup initiated" }</pre> Download

Code

Details

Response headers

```
access-control-allow-credentials: true
content-length: 30
content-type: application/json
date: Fri, 26 Sep 2025 19:58:38 GMT
server: uvicorn
```

Responses

Code	Description	Links
200	Successful Response	No links

Media type

application/json

▼

Controls Accept header.

Example Value

Schema

```
"string"
```

GET

/api/v1/system/audit

Get Audit Logs

Get audit logs

Parameters

Try it out

Name	Description
skip integer (query)	Default value : 0 0
limit integer (query)	Default value : 100 100

Responses

Curl

```
curl -X 'GET' \
'http://127.0.0.1:8000/api/v1/system/audit?skip=0&limit=100' \
-H 'accept: application/json' \
-H 'Authorization: Bearer eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJzdWIiOiIwNDM1YzhINS1iYzI3LTRm
```

Request URL

http://127.0.0.1:8000/api/v1/system/audit?skip=0&limit=100

Server response

Code

Details

200

Response body

```
{
  "user_id": "0435c8e5-bc27-4fff-ac67-c8f538e14936",
  "action": "CREATE",
  "resource_type": "patient",
  "resource_id": "2707a3aa-c798-42af-a94e-40e02ed18a1d",
  "success": true,
  "ip_address": null
},
{
  "id": "583dc86c-5193-4d60-8b33-32697bdb85b7",
  "timestamp": "2025-09-26T19:28:17.562948",
  "user_id": "0435c8e5-bc27-4fff-ac67-c8f538e14936",
  "action": "DELETE",
  "resource_type": "patient",
  "resource_id": "fdd96bb7-08b2-49fb-ad7e-a0f566de59a6",
  "success": true,
  "ip_address": null
},
{
  "id": "86888b96-10e5-47ba-9012-48a014c4cc10",
  "timestamp": "2025-09-26T19:25:58.563040",
  "user_id": "0435c8e5-bc27-4fff-ac67-c8f538e14936",
  "action": "CREATE",
  "resource_type": "patient",
  "resource_id": "166acff6-9d2b-4a50-9be3-da4d241ebf93",
  "success": true,
  "ip_address": null
}
```

Download

Response headers

```
content-length: 1065
content-type: application/json
date: Fri, 26 Sep 2025 19:58:58 GMT
server: uvicorn
```

Responses

Code	Description	Links
200	Successful Response	No links
	Media type	
	application/json	

Code	Description	Links
	Controls Accept header. Example Value Schema <pre>"string"</pre>	
422	Validation Error Media type application/json Example Value Schema <pre>{ "detail": [{ "loc": ["string",], "msg": "string", "type": "string" }]}</pre>	No links

default

^

GET / Root

Root endpoint - API information

Parameters

Try it out

No parameters

Responses

Curl

```
curl -X 'GET' \
'http://127.0.0.1:8000/' \
-H 'accept: application/json'
```



Request URL

http://127.0.0.1:8000/

Server response

CodeDetails

200

Response body

```
{
  "name": "PiroGon Device API",
  "version": "1.0.0",
  "status": "operational",
  "timestamp": "2025-09-26T20:00:11.369084",
  "endpoints": {
    "health": "/api/v1/system/health",
    "docs": "/docs",
    "redoc": "/redoc"
  }
}
```

Download

Response headers

```
content-length: 190
content-type: application/json
date: Fri,26 Sep 2025 20:00:10 GMT
server: uvicorn
```

Responses

Code	Description	Links
200	Successful Response	No links

Media type

application/json

Controls Accept header.

Example Value | Schema

```
"string"
```

GET

/health Health Check



Simple health check endpoint

Parameters

Try it out

No parameters

Responses

Curl

```
curl -X 'GET' \
  'http://127.0.0.1:8000/health' \
  -H 'accept: application/json'
```

Request URL

```
http://127.0.0.1:8000/health
```

Server response

Code	Details
200	Response body <pre>{ "status": "healthy", "timestamp": "2025-09-26T20:00:33.211741" }</pre> Response headers <pre>content-length: 61 content-type: application/json date: Fri,26 Sep 2025 20:00:32 GMT server: uvicorn</pre>

Responses

Code	Description	Links
200	Successful Response Media type application/json Controls Accept header. Example Value Schema <pre>"string"</pre>	No links

Schemas

**Body_complete_session_api_v1_sessions__session_id__complete_put ^**Collapse all **object**

```
test_passed* boolean
stimulation_status* string
notes ^ Collapse all (string | null)
  Any of ^ Collapse all (string | null)
    #0 string
    #1 null
```

Body_create_user_api_v1_auth_create_user_post ^ Collapse all **object**

```
email* string
username* string
full_name* string
password* string
role ^ Collapse all string
  Enum ^ Collapse all array
    #0="admin"
    #1="operator"
    #2="viewer"
  Default="operator"
```

Body_upload_signal_data_api_v1_sessions__session_id__upload__signal_type__Collapse all **object**

```
file* string binary
```

DeviceResponse ^ Collapse all **object**

```
id* string
serial_number* string
name* string
model* ^ Collapse all (string | null)
  Any of ^ Collapse all (string | null)
    #0 string
```

```
#1 null

firmware_version* ^ Collapse all (string | null)

Any of ^ Collapse all (string | null)

#0 string
#1 null

software_version* ^ Collapse all (string | null)

Any of ^ Collapse all (string | null)

#0 string
#1 null

status* ^ Collapse all string

Enum ^ Collapse all array

#0="active"
#1="maintenance"
#2="retired"
#3="offline"

last_calibration* ^ Collapse all (string | null)

Any of ^ Collapse all (string | null)

#0 string date-time
#1 null

next_maintenance* ^ Collapse all (string | null)

Any of ^ Collapse all (string | null)

#0 string date-time
#1 null

registered_at* string date-time
```

DeviceStatus ^ Collapse all string

Enum ^ Collapse all array

```
#0="active"
#1="maintenance"
#2="retired"
#3="offline"
```

HTTPValidationError ^ Collapse all object

detail ^ Collapse all array<object>

Items ^ Collapse all object

loc* ^ Collapse all array<(string | integer)>

Items ^ Collapse all (string | integer)

Any of ^ Collapse all (string | integer)

#0 string

#1 integer

msg* string

type* string

KioskPinRequest ^ Collapse all object

pin* string [4, 6] characters

LoginRequest ^ Collapse all object

email* string email

password* string

OperatorType ^ Collapse all string

Enum ^ Collapse all array

#0="login"

#1="kiosk"

PatientCreate ^ Collapse all object

patient_number* string [1, 50] characters

name* string [2, 100] characters

date_of_birth* ^ Collapse all (string | null)

Any of ^ Collapse all (string | null)

#0 string date-time

#1 null

sex* ^ Collapse all (string | null)

Any of ^ Collapse all (string | null)

#0 PatientSex ^ Collapse all string

Enum ^ Collapse all array

#0="male"

#1="female"

#2="other"

#1 null

medical_record_number ^ Collapse all (string | null)

Any of ^ Collapse all (string | null)

#0 string ≤ 50 characters

#1 null

PatientResponse ^ Collapse all object

patient_number* string [1, 50] characters
 name* string [2, 100] characters
 date_of_birth* > Expand all (string | null)
 sex* > Expand all (string | null)
 medical_record_number > Expand all (string | null)
 id* string
 created_at* string date-time
 updated_at* string date-time

PatientSex ^ Collapse all string

Enum > Expand all array

SessionCreate ^ Collapse all object

session_number* integer [1, 9999]
 threshold_level* number [0, 10]
 tolerance_level* number [0, 10]
 stimulation_level* number [0, 10]
 frequency_level* number [1, 1000]
 pulse_width_level* number [1, 1000]
 operator_name* string [2, 100] characters
 notes ^ Collapse all (string | null)
 Any of ^ Collapse all (string | null)
 #0 string
 #1 null
 patient_id* string matches ^[a-f0-9\-\-]{36}\$
 device_id* string matches ^[a-f0-9\-\-]{36}\$
 operator_type ^ Collapse all string
 Enum ^ Collapse all array
 #0="login"
 #1="kiosk"
 Default="login"

SessionResponse ^ Collapse all object

```
session_number* integer [1, 9999]
threshold_level* number [0, 10]
tolerance_level* number [0, 10]
stimulation_level* number [0, 10]
frequency_level* number [1, 1000]
pulse_width_level* number [1, 1000]
operator_name* string [2, 100] characters
notes ^ Collapse all (string | null)
  Any of ^ Collapse all (string | null)
    #0 string
    #1 null
id* string
start_time* string date-time
end_time* ^ Collapse all (string | null)
  Any of > Expand all (string | null)
status* ^ Collapse all string
  Enum ^ Collapse all array
    #0="active"
    #1="completed"
    #2="aborted"
    #3="error"
patient_id* string
device_id* string
operator_type* ^ Collapse all string
  Enum ^ Collapse all array
    #0="login"
    #1="kiosk"
test_passed* ^ Collapse all (boolean | null)
  Any of ^ Collapse all (boolean | null)
    #0 boolean
    #1 null
stimulation_status* ^ Collapse all (string | null)
  Any of ^ Collapse all (string | null)
    #0 string
    #1 null
```

SessionStatus ^ Collapse all **string****Enum** ^ Collapse all **array**

#0="active"

#1="completed"

#2="aborted"

#3="error"

TokenResponse ^ Collapse all **object****access_token*** **string****token_type** ^ Collapse all **string****Default**="bearer"**expires_in*** **integer****user*** **object****UserRole** ^ Collapse all **string****Enum** ^ Collapse all **array**

#0="admin"

#1="operator"

#2="viewer"

ValidationError ^ Collapse all **object****loc*** ^ Collapse all **array<(string | integer)>****Items** ^ Collapse all **(string | integer)****Any of** ^ Collapse all **(string | integer)**#0 **string**#1 **integer****msg*** **string****type*** **string**